

Facial Recognition

Author

Name: Aqsa mahmood Email: aqsa.mahmood@city.ac.uk

Hardware Specification

This project was run on the following hardware specification: | Hardware | Description
| | ----- | ----- | | System Model: | MS-7B93 | | System Type | x64-based PC | |
Processor(s): | AMD64 Family 23 Model 113 Stepping 0 AuthenticAMD ~3800 Mhz | |
Graphics | NVIDIA GeForce RTX 2070 SUPER | | BIOS Version: | American Megatrends Inc.
1.00, 14/06/2019 | |Total Physical Memory: | 16,332 MB|

Code Structure

```
GANKneeMRI2
|  README.md
|  environment.yml.yml
|  setup.bat
|  setup.sh
└── src
    |  └── models
    |  |  mrnet.py
    |  |  resnet_G_and_C.py
    |  └── utils
    |  |  data_load.py
    |  |  fid.py
    |  |  knee_mri_viewer.ipynb
    |  |  logger.py
    |  |  utilities.py
    |  └── training
    |  |  baseline_model_training.py
    |  |  wasserstein_GP_256.py
    |  └── traininer
    |  |  basenet_model_task.py
    |  |  wgan_task.py
    |  └── bash_scripts
    |  |  base-mode-train-cloud.bat
    |  |  base-mode-train-cloud.sh
    |  |  wgan-model-train.bat
    |  |  wgan-model-train.sh
    |  |  mrnet-model-train.bat
    |  |  mrnet-model-train.sh
    |  └── model_weights
    |  |  ....
    |  |  ....
    |  |  Dockerfile
    |  |  Evaluating_classificationPerf .py
    |  |  mrNet_main.py
    |  |  setup.py
```

```
| | testing_GAN_model.py
| | wgan_main.py
```

Python File/Folder Description

- testing_GAN_model.- this file essentially allows test mri scans to be generated from the trained WGAN models
- mrNet_main.py - this file essentially allows the training of the baseline model
- wgan_main.py - this file essentially allows the training of WGAN
- utilities.py - this file contains general utility functions used throughout the files
- dataset_loader.py - contains functions to allow loading the MRI Dataset
- mrnet.py - contains the baseline model
- resnet_G_and_C.py - contains WGAN Models Discriminator and Generator
- setup.py - this file is for Google Cloud purposes and is not yet functional. The contents of this can be ignored.
- Dockerfile - this file is for Google Cloud purposes and is not yet functional
- /trainer - this directory is for Google Cloud purposes and is not yet functional

Model Downloads

The 'model_weights' directory located in src folder contains weights of the trained WGAN models

Installation

First, install [Python 3.7](#) and [Anconda](#). Then we need to setup the anaconda environment. On a linux environment, navigate to the FacialRecognition project, run the following on command line as super user:

```
$ ./setup.sh
```

This script should create the conda environment, with all development dependencies, and activate the environment 'project2020' | Dependency | Version | | ----- | ----- | | python | 3.7.6 | | pytorch | 1.2.0 | | pillow | 6.1 | | torchvision | 0.4.0 | | fastai | 1.0.57 | | ipykernel: | 4.6.1 | | pytest | 3.6.4 | | bqplot | | scikit-learn | 0.19.1 | | pip | 19.0.3 | | cython | 0.29.1 v | papermill | 1.2.0 | | black | 18.6b4 v | ipywebtrc | | lxml | 4.3.2 | | pre-commit | 1.14.4 | | pyyaml | 5.1.2 | | requests | 2.22.0 v | einops | 0.1.0 | | cytoolz | | seaborn | | mtcnn | | scikit-image | | decord | 0.3.5 | | nvidia-ml-py3 | | nteract-scrapbook | | azureml-sdk[notebooks,contrib] | 1.0.30 | | facenet-pytorch | | opencv-python | 3.4.2.16 | | opencv-contrib-python | 3.4.2.16 | | jupyter | 1.0.0 | | spyder | 4.1.3 |

On a windows environment, can run the following script via Anaconda Prompt :

```
$ setup.bat
```

Note that setting up the environment does take some time. Once, the the script has finished, the environment 'project2020'. On the same Anaconda Prompt, spyder can be launched via the following command:

```
$ spyder
```

Running Application

The python files of interest is **wgan_main.py** and **mrNet_main.py**. These python files can be invoked via scripts located in the 'bash_scripts' directory. The wgan-model-train.bat will run the wgan_main.py file, and will start WGAN training. The batch file can be invoked via Anaconda Prompt with the following example command:

```
$ wgan-model-train.bat --base-results-dir="D:/results" --data-dir="C:/Users/adds0/OneDrive - University of Glasgow/Documnets/courses/Msc/Project/MRI/MRNet/resources" --plane=axial --diagnoses=ACL --dataset-type=train
```

The above argument call would train a WGAN model for the dataset that is of axial plane and which have a diagnosis of ACL. Note that the first argument '--base-results-dir' is the results directory path, and the second argument, '--data-dir', is the location of the original dataset.

The wgan-model-train.sh can be run from Anaconda Prompt by invoking the following command:

```
$ wgan-model-train.sh --base-results-dir="D:/results" --data-dir="C:/Users/adds0/OneDrive - University of Glasgow/Documnets/courses/Msc/Project/MRI/MRNet/resources" --plane='axial' --diagnoses='ACL' --dataset-type='train'
```

***Note before invoking the wgan-model-train.sh, it needs to be made executable with the following command:**

```
$ chmod u+x wgan-model-train.sh
```

If one wants to train a WGAN model for a different plane, then one can specify 'sagittal' or 'coronal'. Alternatively, if WGAN need to be trained for a different diagnoses, the following options are available:

- 'Abnormal' - only selects scans that are only abnormal
- 'Meniscus' - only selects scans that are only Meniscus
- 'All' - selects scans that are all ACL, Meniscus and Abnormal
- 'None' - selects scans that are normal

The mrnet-model-train.bat will run the mrNet_main.py file, and will start base model training with or without the augmented dataset. The batch file can be invoked via Anaconda Prompt with the following example command:

```
$ mrnet-model-train.bat --base-results-dir="D:/results" --data-dir="C:/Users/adds0/OneDrive - University of Glasgow/Documnets/courses/Msc/Project/MRI/MRNet/resources" --plane=axial --dataset-type=train --epochs=100
```

The above argument call would train a Base MRNet model for the dataset that is of axial plane without an augmented dataset. Note that the first argument '--base-results-dir' is the results directory path, and the second argument, '--data-dir', is the location of the original dataset. To include training with an augmented dataset, specify '--augmented-data-dir' with the directory which contains the augmented dataset. Ensure the directory structure is similar to the contents within '--data-dir'. Again, if one wants to train the baseline model for a different plane, then one can specify 'sagittal' or 'coronal'.